

**PRODUCT**  
**CATALOG**  
2017 / 2018

3,75


Qp

ENG  
v3.0



Implant Reliability Quality Warranty Surgical tool Surgical kits Prosthetic About Drills Surgical prosthetic

Long Skills Education Plastic sleeve Wide range of sizes Anatomic SGS De Bone type Ratchet Hex Driver Favorable prices Conus platform Ball attachment Angular Tilling Cap Hex Titanium Zirconium Knowledge Surgical tools



## TABLE OF CONTENTS



Company profile .....	4 - 5
Surface treatment .....	6 - 11
Packaging system.....	12
Sign index .....	13
Implant systems.....	14 - 15
Surgical drills .....	26-27
Surgical kits.....	28 - 31
Stopper drills kits.....	32 - 33



P7S dental implant .....	22 - 23
P9S dental implant.....	24
Surgical tools for P7S/P9S .....	25

Surgical Tools .....	17
P7 dental implant.....	18 - 19
P1 dental implant.....	20 - 21
Healing abutments .....	36
Impression transfers .....	37
Analog .....	38
Straight abutments.....	39
Angular abutments.....	40
Anatomic angular abutments.....	41
Zirconium abutments .....	42 - 43
Temporary abutments.....	44
Abutments for casting.....	45
CAD-CAM .....	46 - 47
Overdenture ball attachments .....	48 - 49
Abutments for immediate loading .....	50 - 51
Overdenture S-Lock abutments.....	52 - 53
Overdenture Easy-Fix abutments.....	54 - 55
Overdenture Smart - Lock abutments.....	56 - 57
Straight multi-unit abutments .....	58 - 59
Multi-base abutments .....	60 - 61
The-One multi-unit abutments.....	62 - 63



## Company profile

**S**GS International Ltd. had been incorporated in Liechtenstein, Schaan with aim to offer the customers a comprehensive range of innovative, science-based dental solutions. SGS Dental Implant System is not only a dental product, but it is about art and experience in clinical practise as well.

Thanks to qualitative Swiss technology SGS Dental Implant System provides excellent products on competitive prices with strong support. We use the most precise tools and equipments with the best raw material, that is why SGS became one of market leaders in the dental implant industry and nowadays our system is used by thousands of professionals all around the world.

All our implants are made of biocompatible material - medical grade 5 titanium alloy with extraordinary coating. Each our product meet the most stringent international quality standards, inspected with very precise triple quality control. SGS International Ltd. has obtained the CE by the European Notified Body CE 1979, European Directive 93/42/EEC - Annex II, Section 4 and are also certified with the quality standard EN ISO 9001:2008 and ISO 13485:2012 on devices dental-implant. SGS International Ltd. is also registered by NQA quality management system. SGS Medical Devices have been cleared for marketing at the US market as well.

2009 is a milestone in our company's life opening a branch

and afterwards a brand new 2500 sqm site with new sterile packaging, logistics and worldwide distribution centre in European Union - Hungary, Budapest. Situated at the very heart of Europe, having hence excellent conditions for expanding, nowadays we have been providing service for our customers in more than 20 countries in the world. Besides of high quality European products SGS Dental Implant System gives best quality service, guarantees affordable price and takes care of all Partners, always intending to be in close business relationship with them. One of our strengths is the worldwide Distribution Network. We are very proud of our well educated distributors providing high-level professional support to their partners.

Hold in high esteem practical knowledge and education in 2013 we have opened our SGS Medical Centre, where medical educational trainings, implantology courses can be provided for oral surgeons visiting us from different part of the world. Thanks to these trainings, courses doctors have well-founded knowledge and experience in SGS Dental Implant System. They are able to work with qualitative and innovative European dental implant system, moreover company SGS guarantees maximum support, help and fast service.

We have been developing and trust to achieve measureless success with all the desired goals with our partners, as we have

**ALL SOLUTIONS FOR A PERFECT SMILE.**





## Titanium material

All SGS Dental Implants are made from European "bio friendly" medical titanium alloy, grade 5 - Titanium 6AL-4V, (signifying the Titanium alloy containing 6% Aluminium and 4% Vanadium alloy). Titanium's special property of fusing to bone, called osseointegration ("osseo" – bone; "integration" – fusion or joining with), is the biological basis of dental implant success.

When teeth are lost, the bone that supported those teeth is lost too, so placing dental implants stabilizes bone, preventing its loss. Along with replacing lost teeth, implants help maintain the jawbone's shape and density. Dental implants help you eat, chew, smile, talk and look completely natural.



## SGS sterilization procedure

SGS has a strict sterilisation process, undergoing a continuous laboratory control and report manual management. All our implant undergo a clean-room packaging channel and then introduced to gamma radiation. We keep a hard control on the steps, our colleagues practise and each dosage of implants is randomly inspected in SGS laboratories.



## Product guarantee

SGS Dental Implant takes responsibility for all of its implants and gives life-time warranty when used in accordance with the supplied instructions for use and the company protocols.



Our implants  
available in

**28**  
countries



*Out of 100 implantations less than two have negative result because of our special smart surface treatment, unique implant design and our precious research.*

## SGS Dental Implants



Dental implants are an effective, safe and excellent solution to the problems resulting from missing teeth, look and feel like your own teeth. That is why in the last 25 years dental implants have changed the face of dentistry.

A dental implant is actually a replacement for the root of tooth, secured in the jawbone and not visible once surgically placed. It is used to secure crowns (the parts of teeth seen in the mouth), bridgework or dentures by a variety of

means. Nowadays good and safe implant is not only a prosthetic, but aesthetic question as well. Extremely important that our implants are made of lightweight, strong and biocompatible titanium. Biocompatible word means not rejected by the body, furthermore we have excellent surface treatment which reduces healing time.



ISO 9001:2008  
ISO 13485:2012

SGS's products are cleared for marketing in the USA.\*

*\*Some products may not be available in the USA.*




## A new microstructured bioactive antibacterial surface for implants!

SGS Dental implements famous SBTC® coating for its dental implants:

SBTC® is a known worldwide type of dental implant coating, having outstanding performances in dental implantation practice.


### Advantages of the SBTC® coating

- + Faster and better healing
- + Complex surface design with significant surface enlargement
- + High hydrophilic reaction with blood
- + Increased primary stability with reduced healing time
- + Active support of bone attachment
- + Higher application security
- + Possible diversification of indications (early loading/immediately loading)
- + Prevention of spontaneous oxidation of the titanium surface through CaP-coating
- + Higher osteoconductivity of the surface
- + Outstanding biocompatibility
- + Thin coating
- + Microcrystalline structure, large open surface
- + High solubility and controlled resorption area
- + Complete coverage of porous surfaces and complex implant geometries
- + Microporosity with high capillary effect on body fluids






SEM magnificationX25: clean uniform surface with no contaminations




SEM magnificationX2000: clear uniform crystalline structure of SBTC® type




EDS spectrum: Calcium-Phosphorus-Oxygen presence adequate to chemical composition of SBTC® coating

### SBTC coated dental implants of SGS Dental have all substantial features of the SBTC type Ca-P coating:


- + Clean uncontaminated uniform surface
- + Unique SBTC crystalline structure of Ca-P brushite particles' coating on dental implants surface
- + Ca-P-O chemical composition of the coating layer approving its SBTC origin



Osteocalcin, 14 days



Osteocalcin, 30 days




### Summary

Implants with SBTC surface showed a significant increase in Bone-Implant-Contact (BIC) in the spongiosa area between 14 and 30 days. In the further course within the SBTC remodeling BIC-Data in the range between 40-60% arises, which conforms to the data described in the literature. After 30 days the osteocalcin –expression too was significant increased by the implants with the SBTC-surface.

### Description of the biological properties of the coating


The SBTC® coating is a bioactive calcium phosphate coating that supports the adhesion of osteoblast cells and simultaneously promotes their proliferation. The cells demonstrate good adhesion and a typical morphology for osteoblast. Under the scanning electron microscope the integration of the cells into the material is clearly visible.



Bone tissue formation on SBTC®



Human osteoblasts on SBTC®




Osteoblast MG 63 cells on SBTC®

The SBTC® coating consists of two calcium phosphate phases with different solubilities. The more easily soluble outer calcium phosphate phase, brushite, occurs in natural bone as an intermediate stage during calcification of new bone tissue. When brushite dissolves, calcium and phosphate ions are released in a high concentration, and they are the cause of fast contact osteogenesis and the high mineralization rate. Brushite is therefore in a position to stimulate the body to its own bone synthesis in the short term, and to accelerate the osseointegration of the implants, particularly in the primary phase. The inner calcium phosphate phase, the fine crystalline hydroxyapatite, is resorbed more slowly and releases ions that promote the formation of new bone over a longer period. The SBTC® coating is fully resorbed over a period of 6-12 weeks after implant placement and is simultaneously replaced by newly formed bone tissue, with the ultimate result that an optimum bond between bone and implant has been formed in place of the coating. This osteoinductive property combined with the controlled resorption is the primary advantage of the bioactive SBTC® coating.

## Differentiation of cells in vitro under the influence of SBTC®


The influence of the SBTC® coating on cell differentiation was examined by a Co culture of the hFOB1.19 osteoblast cell line with TPS/SBTC®-coated platelets of TiAl6V4. The osteoblast-specific collagen synthesis was analyzed at various points during incubation. The result after 6 days and after 10 days of incubation showed increased collagen synthesis on the SBTC®-coated test bodies.



## Mineralization in vitro under the influence of SBTC®

The influence of the SBTC® coating in the mineralization was analyzed by incubating test bodies coated with SBTC® in cell culture medium (DMEM) for seven days. The extract was added to a confluent cell layer and the mineralization was confirmed by van Kossa staining. With

van Kossa staining mineralized areas are stained black. Figure 10 shows the difference between the control medium and the SBTC® extract. While a slight mineralization could be confirmed in the cells in the control medium, strong mineralization could be confirmed with the SBTC® extract. This indicates that the calcium-phosphate phases in the SBTC® coating stimulate the mineralization of human osteoblasts.




## Immunological reactivity in vitro under the influence of SBTC®

The effect of the SBTC® coating on the immunological reactivity was analyzed by the release of interleukin 1 $\beta$  (IL-1 $\beta$ ). IL-1 $\beta$  is a typical enzyme, which is released during the early inflammation phase and influences bone resorption. The test was conducted with monocytes and macrophages of the mouse cell line J-774A.1, which were cultured either with control bodies (TiAl6V4/TPS) or with SBTC®-coated test bodies. After three days of culture the expression of IL-1 $\beta$  was analyzed. The SBTC®-coated samples in comparison with the uncoated control demonstrated a significant reduction of IL-1 $\beta$  release. This means that the coating with SBTC® triggers virtually no inflammation parameters, and as a result can be classified as very biocompatible.






### Immunological Reactivity



### Analysis of the protein adsorption in vitro under the influence of SBTC®

The protein adsorption or immobilization of proteins at the implant surface is, clinically viewed, an important step in the osseointegration of implants. To determine the protein adsorption SBTC®-coated test bodies and uncoated control bodies were incubated in fetal calf serum for several hours. After various incubation times (1h and 4h) the protein adsorption on the different test bodies was analyzed. As can be seen in the next diagram, the coating with SBTC® significantly increased protein adsorption in comparison with the uncoated test bodies.



after 8h


after 30h

after 48h

### Conclusion

The electrochemical deposition of the SBTC® coating ensures complete coverage of porous implant surfaces and complex implant shapes. Unlike the highly crystalline, poorly soluble HA coatings in the plasma-spray process, the electrolytic coating technology yields a fine crystalline structure. The process eliminates hard particles and delamination of the coating. The almost vertical adjacent calcium phosphate crystals and the associated large open surface give the implant surface a high capillary effect on blood and ensure adsorption and immobilization of relevant growth factors. The controlled process of coating degradation correlates with simultaneous formation of new bone, which takes place immediately on the porous

### Protein Adsorption



### In vitro precipitation tests with SBTC®

In in vitro tests SBTC®-coated test bodies were colonized with osteoblast cells of the cell line MG-63 and cultured in cell culture medium for 48 hours. After 30 hours in the culture a fine crystalline precipitation could be observed on the coating surface. The cells on the SBTC® surface were partially covered with the precipitate. As could be shown by EDX analyses, the precipitate was also a calcium-phosphate compound. After 48 hours the cells were completely covered. Visualization of the actin cytoskeleton of the bone cells showed that the morphology of the cells remained virtually unchanged during the reprecipitation.


Comparison test bodies that were coated with hydroxyapatite only did not show this behavior. The in vitro results lead to the conclusion that there is a precipitation on the surface during dissolution of the coating, in particular because of the presence of the easily dissolved brushite phase. It can be concluded that these processes also take place in the body and there is therefore a calcium phosphate phase directly on the SBTC® surface in the body.

implant surface. This results in increased bone deposition and the option of early mechanical loading. This gives the SBTC® coating improved osseointegration and it can be considered a further development of the plasma-sprayed HA coatings, which retains the good bioactive properties and eliminates the remaining potential long-term risks. The unique two-layer design of the SBTC® coating is perfectly adapted to the healing process of the bone. The SBTC® coating is resorbed within 6-12 weeks by a controlled mechanism and completely replaced by new bone tissue. This means that the coating only remains in place until the implant is healed in position and has formed a fixed connection between the implant surface and the surrounding bone tissue.




## Animal-experimental tests

Animal tests were conducted on the minipig animal model at the laboratory for biomechanics and experimental orthopedics of the University of Mannheim to investigate the formation of new bone with SBTC®-coated implants. SBTC®-coated titanium pins were press-fitted into 21 animals. The follow-up period was 12 weeks. The results of the study show significantly increased bone deposition with the SBTC®-coated implants and subsequently significantly better anchorage of the implant in the early postoperative phase.



Shear strength of different surfaces after implant placement




Bone deposition of different surfaces after implant placement

## Effectiveness of the SBTC® coating in the animal model

In this animal experiment the osseointegration of test implants with the TPS surface was compared to

implants with the TPS/SBTC® coating. The implants were placed in the maxilla of domestic pigs (*Sus scrofa domestica*). The direct bone contact in both test groups was analyzed six weeks after implant placement. The results of the test show significant differences in the bone contact between the two groups. The average bone contact for the control implants was 49.8%, while a direct bone contact of 73% was measured for the SBTC®-coated implants.



Bone contact values of the different surfaces

A high proportion of bone was found, particularly between the threads of implant, and a clear osteoconductive effect as a result of the presence of the SBTC® coating. The SBTC® coating was almost completely resorbed during the study period of six weeks.

No reactions to foreign bodies were detected, which is another indication of the very good biocompatibility of the surface.


## Study of SBTC®-coated implants in the canine model

The goal of the study was to determine the effect of the SBTC® coating on osseointegration over extended periods with immediate loading. Implants with different surfaces were placed in the mandibles of dogs (beagles). The surfaces of the implants consisted of a TPS surface, a plasma-sprayed HA surface and a TPS+SBTC® surface. The implants were immediately restored with a crown and placed under immediate loading. The crowns were not in contact with neighboring teeth or other implants. The follow-up period was seven months. The results of the trial show that the SBTC® coating was fully resorbed after seven months and had been replaced by newly formed bone tissue.



In contrast, fragmentation of the coating and unhomogenous resorption could be observed with the HA-plasma-sprayed surface. Isolated HA particles were also found.

The SBTC®-coated implants also demonstrated the highest bone deposition density. However, the difference between the surfaces decreased with increasing implant placement time



flakes, while the electrochemically deposited CaP coating had rod-shaped crystals with a hexagonal cross-section.

The histological analyses after six weeks showed bone growth along the surfaces. On the electrochemically deposited CaP coating the significantly largest BIC values were measured and compared with the rough and biomimetically deposited CaP surfaces. The study showed that the electrochemically deposited CaP coating appears to improve osseointegration, and as a result can ensure a long-term and stable fixation of the implants in the bone tissue.

## Clinical results

A multicentric study, which included universities and private practices, investigated PITT-EASY implants (Oraltronic) with SBTC® coating (FBR surface on a porous TPS surface). The implants were placed in the maxilla and the mandible. The study protocol included immediate loading. A total of 156 implants were placed in 62 patients, with 40 implants placed in fresh extraction alveoli. After 6 months 8 implants in 6 patients had been lost, 6 in the mandible and 2 the maxilla. After 6 months under load 98% of the implants were osseointegrated and functional.



## Effect of differently applied CaP coatings on the osseointegration of implants

The effect of differently applied CaP coatings on the osseointegration of titanium implants was investigated in the animal model. The study included three groups with different surface modifications. Group 1 had a rough surface, group 2 had a biomimetic CaP coating and group 3 had an electrochemically deposited CaP coating. A total of 36 implants were placed in the tibias of 18 rabbits. The study period was 6 and 12 weeks. The influence of the different implant surfaces on the osseointegration was analyzed. REM images of the different surfaces were prepared and analyzed. On the biomimetically deposited CaP coating the crystals were arranged as




# SGS DENTAL PACKAGING SYSTEM



Remove the implant from the plastic holder:









After you removed the implant from the plastic holder:



Surgical tool mount




## Sign Index

<b>REF</b>	Catalogue number
<b>LOT</b>	Batch number
	Manufacturer
	Use by date
	Do not use if package is opened or damaged
	Do not re-use
	Do not resterilize
	Attention, see instruction for use
<b>RxOnly</b>	Prescription use only
<b>STERILE R</b>	Sterilized
<b>CE</b>	Conformity marking


Color codes:



With color codes it is easy to identify the implant's diameter.







3,75

**PLATFORM  
3.75**

**ONE-PIECE**  
IMPLANT

**ONE-PIECE**






## P1

Screw Type Groovy Implant

Short pitch, slightly aggressive implant, which is specially for hard bone (mandibula), with platform shift and internal hex.

Available sizes:  
Ø 3.2, Ø 3.75, Ø 4.2, Ø 5, Ø 6



## P7

Conical Groovy Implant

High pitch, slightly aggressive implant, which is suitable for immediate loading. It is characterized by self tapping screw with platform shifting and internal hex.

Available sizes:  
Ø 3.2, Ø 3.75, Ø 4.2, Ø 4.5, Ø 5, Ø 6



## P7S

Integral Groovy Implant

One-piece implant. It has the same parameters like P7. The abutment is integrated.

Available sizes:  
Ø 3.0, Ø 3.2, Ø 3.75, Ø 4.2, Ø 5, Ø 6



## P9S

Thin Integral Implant


One-piece esthetic implant.

Available sizes:  
Ø 2.4, Ø 3.0, Ø 3.2








## PLATFORM 3.75

---



Please note the image shown is for illustration purposes only!









## Instrument tools and drivers

					
Product code	K5 - 0 mm	K5 - 9 mm	K5 - 15 mm	K8 - 23 mm	K8 - 28 mm
Ref. number	B50	B59	B515	B823	B828
Length	0 mm	9 mm	15 mm	23 mm	28 mm
Description	Driver for implants			Adaptor for implants	
Material	Stainless steel				

## Prosthetic tools for the abutments

							
Product code	K1 - 9 mm	K1 - 15 mm	K2 - 9 mm	K2 - 15 mm	K9 - 23 mm	DG	K7
Ref. number	B19	B115	B29	B215	B923	C125	B7
Length	9 mm	15 mm	9 mm	15 mm	23 mm		18 mm
Description	Driver for abutments		Hand driver for abutments		Adaptor for abutments	Depth checker	Retrieval screw
Material	Stainless steel					Titanium	

## Implantology tools

								
Product code	K3D/K3M	D9	D17	R8	R8T	M+	DP	K10
Ref. number	B30001/B30002	C87	C88	C97	C98	C96	C8920/C8927	B0
Description	Adaptor for driver	Handle	Tissue punch	Ratchet wrench	Ratchet torque	Depth gauge	Parallel pin	Technical driver
Material	Stainless steel							

**3,2**

Length

8 mm

10 mm

11,5 mm

13 mm

16 mm

Ref. number

M02328

M023210

M023211

M023213

M023216

**3,75**

Length

8 mm

10 mm

11,5 mm

13 mm

16 mm

Ref. number

M02378

M023710

M023711

M023713

M023716

**4,2**

Length

6 mm

8 mm

10 mm

11,5 mm

13 mm

16 mm

Ref. number

M02426


M02428

M024210

M024211

M024213

M024216



- Marker drill - to be used to make only a mark
- Throughout entire implant's length
- Drill only through the cortical bone, should not be used to full depth.  
If the cortical bone is hard (D1), you may use this drill as a countersink.

An additional 0,8 - 1,0 mm must be added to the length of the drill to account for the angled cutting up.  
Procedure recommended by SGS cannot replace the judgment and the experience of the surgeon!

**4,5**

Length

6 mm

8 mm

10 mm

11,5 mm

13 mm

16 mm

Ref. number

M02456

M02458

M024510

M024511

M024513

M024516

**5**

Length

6 mm

8 mm

10 mm

11,5 mm

13 mm

16 mm

Ref. number

M0256

M0258

M02510

M02511

M02513

M02516

**6**

Length

6 mm

8 mm

10 mm

11,5 mm

13 mm

16 mm

Ref. number

M0266


M0268

M02610

M02611

M02613

M02616



- Marker drill - to be used to make only a mark
- Throughout entire implant's length
- Drill only through the cortical bone, should not be used to full depth.  
If the cortical bone is hard (D1), you may use this drill as a countersink.

An additional 0,8 - 1,0 mm must be added to the length of the drill to account for the angled cutting up.  
Procedure recommended by SGS cannot replace the judgment and the experience of the surgeon!



**3,2**


Length	8 mm	10 mm	11,5 mm	13 mm	16 mm
Ref. number	M01328	M013210	M013211	M013213	M013216

**3,75**

Length	8 mm	10 mm	11,5 mm	13 mm	16 mm
Ref. number	M01378	M013710	M013711	M013713	M013716

**4,2**

Length	8 mm	10 mm	11,5 mm	13 mm	16 mm
Ref. number	M01428	M014210	M014211	M014213	M014216








- Marker drill - to be used to make only a mark
- Throughout entire implant's length
- Drill only through the cortical bone, should not be used to full depth.  
If the cortical bone is hard (D1), you may use this drill as a countersink.


An additional 0,8 - 1,0 mm must be added to the length of the drill to account for the angled cutting up.  
Procedure recommended by SGS cannot replace the judgment and the experience of the surgeon!

**5**

						
Length	6 mm	8 mm	10 mm	11,5 mm	13 mm	16 mm
Ref. number	M0156	M0158	M01510	M01511	M01513	M01516

**6**

						
Length	6 mm	8 mm	10 mm	11,5 mm	13 mm	16 mm
Ref. number	M0166	M0168	M01610	M01611	M01613	M01616



- Marker drill - to be used to make only a mark
- Throughout entire implant's length
- Drill only through the cortical bone, should not be used to full depth.  
If the cortical bone is hard (D1), you may use this drill as a countersink.

An additional 0,8 - 1,0 mm must be added to the length of the drill to account for the angled cutting up.  
Procedure recommended by SGS cannot replace the judgment and the experience of the surgeon!

**3**

P7S  
NARROW

Length	10 mm	11,5 mm	13 mm	16 mm
Ref. number	O05310	O05311	O05313	O05316


**3,2**

P7S  
NARROW

Length	8 mm	10 mm	11,5 mm	13 mm	16 mm
Ref. number	O05328	O053210	O053211	O053213	O053216

**3,75**

Length	8 mm	10 mm	11,5 mm	13 mm	16 mm
Ref. number	O05378	O053710	O053711	O053713	O053716



- Marker drill - to be used to make only a mark
- Throughout entire implant's length
- Drill only through the cortical bone, should not be used to full depth.  
If the cortical bone is hard (D1), you may use this drill as a countersink.

An additional 0,8 - 1,0 mm must be added to the length of the drill to account for the angled cutting up.  
Procedure recommended by SGS cannot replace the judgment and the experience of the surgeon!

**4,2**

					
Length	8 mm	10 mm	11,5 mm	13 mm	16 mm
Ref. number	O05428	O054210	O054211	O054213	O054216

**5**

						
Length	6 mm	8 mm	10 mm	11,5 mm	13 mm	16 mm
Ref. number	O0556	O0558	O05510	O05511	O05513	O05516

**6**

						
Length	6 mm	8 mm	10 mm	11,5 mm	13 mm	16 mm
Ref. number	O0566	O0568	O05610	O05611	O05613	O05616




Drill Speed (RPM)    700-800    700-900    500-700    500-700    400-700    400-600    400-600    400-600    300-500

Diameter    Ø 1.5    Ø 2.0    Ø 2.5    Ø 2.8    Ø 3.2    Ø 3.7    Ø 4.0    Ø 4.5    Ø 5.5


**Ø 4.2 mm**

**Ø 5 mm**

**Ø 6 mm**



-  Marker drill - to be used to make only a mark
-  Throughout entire implant's length
-  Drill only through the cortical bone, should not be used to full depth.  
If the cortical bone is hard (D1), you may use this drill as a countersink.

An additional 0,8 - 1,0 mm must be added to the length of the drill to account for the angled cutting up.  
Procedure recommended by SGS cannot replace the judgment and the experience of the surgeon!

				
Length	10 mm	11,5 mm	13 mm	16 mm
Ref. number	R062410	R062411	R062413	R062416

	<div>3,0</div> <div>P95 NARROW</div>				
Length	10 mm	11,5 mm	13 mm	16 mm	
Ref. number	R06310	R06311	R06313	R06316	

	<div>3,2</div> <div>P95 NARROW</div>				
Length	10 mm	11,5 mm	13 mm	16 mm	
Ref. number	R063210	R063211	R063213	R063216	



- Marker drill - to be used to make only a mark
- Throughout entire implant's length
- Drill only through the cortical bone, should not be used to full depth.  
If the cortical bone is hard (D1), you may use this drill as a countersink.

An additional 0,8 - 1,0 mm must be added to the length of the drill to account for the angled cutting up. Procedure recommended by SGS cannot replace the judgment and the experience of the surgeon!






## Instrument tools and drivers



Product code	K5 - 0 mm	K5 - 9 mm	K6	K8 - 23 mm	K8 - 28 mm
Ref. number	B50	B59	ORB06	B823	B828
Length	0 mm	9 mm	16 mm	23 mm	28 mm
Material	Stainless steel				
Instructions					



## Analog and transfer for P7S & P9S


Product code	A1 - P9S/P7S-N	T1 - P9S/P7S-N	A1 - P7S	T1 - P7S
Ref. number	OR08	OR77	O08	O77
Dimensions:	D: 3 mm	D: 3 mm	D: 3.6 mm	D: 3.6 mm
Material	Stainless steel	Plastic	Stainless steel	Plastic
Instructions	⚠ For P9S and P7S Ø3.0 / Ø3.2		⚠ For P7S Ø3.75 / Ø4.2 / Ø5 / Ø6	






































Straight Drills								
								
Diameter	Ø 2.0	Ø 2.5	Ø 2.8	Ø 3.2	Ø 3.7	Ø 4.0	Ø 4.5	Ø 5.5
Product code	DR - 2.0	DRF - 2.5	DRF - 2.8	DRF - 3.2	DRF - 3.7	DRF - 4.0	DRF - 4.5	DRF - 5.5
Ref. number	C992	C9025	C9028	C9032	C9037	C9040	C9045	C9055
Material	Stainless steel							
Information	Biocompatible diamond link carbon coating							





Conical Drills								
								
Diameter	Ø 2.0	Ø 2.5	Ø 2.8	Ø 3.2	Ø 3.7	Ø 4.0	Ø 4.5	Ø 5.5
Product code	DR - 2.0	DRC - 2.5	DRC - 2.8	DRC - 3.2	DRC - 3.7	DRC - 4.0	DRC - 4.5	DRC - 5.5
Ref. number	C992	C9125	C9128	C9132	C9137	C9140	C9145	C9155
Material	Stainless steel							

Short Drills								
								
Diameter	Ø 2.0	Ø 2.5	Ø 2.8	Ø 3.2	Ø 3.7	Ø 4.0	Ø 4.5	Ø 5.5
Product code	DRS - 2.0	DRS - 2.5	DRS - 2.8	DRS - 3.2	DRS - 3.7	DRS - 4.0	DRS - 4.5	DRS - 5.5
Ref. number	C1002	C10025	C10028	C10032	C10037	C10040	C10045	C10055
Material	Stainless steel							





Stopper Drills							
6 mm							
Diameter	Ø 2.0	Ø 2.5	Ø 2.8	Ø 3.2	Ø 3.7	Ø 4.0	Ø 4.5
Ref. number	C10126	C101256	C101286	C101326	C101376	C101406	C101456
8 mm							
Diameter	Ø 2.0	Ø 2.5	Ø 2.8	Ø 3.2	Ø 3.7	Ø 4.0	Ø 4.5
Ref. number	C10128	C101258	C101288	C101328	C101378	C101408	C101458
10 mm							
Diameter	Ø 2.0	Ø 2.5	Ø 2.8	Ø 3.2	Ø 3.7	Ø 4.0	Ø 4.5
Ref. number	C101210	C1012510	C1012810	C1013210	C1013710	C1014010	C1014510
11.5 mm							
Diameter	Ø 2.0	Ø 2.5	Ø 2.8	Ø 3.2	Ø 3.7	Ø 4.0	Ø 4.5
Ref. number	C101211	C1012511	C1012811	C1013211	C1013711	C1014011	C1014511
13 mm							
Diameter	Ø 2.0	Ø 2.5	Ø 2.8	Ø 3.2	Ø 3.7	Ø 4.0	Ø 4.5
Ref. number	C101213	C1012513	C1012813	C1013213	C1013713	C1014013	C1014513

Other Drills			 <div>TREPHINE DRILL 4 mm</div> <div>D1 D2</div>	 <div>TREPHINE DRILL 5 mm</div>
Product code	DR - 1.5	DR-E	DR-T 4/5	DR-T 5/6
Ref. number	C9215	C93	C9445	C9456
Dimensions			D1: 5 mm   D2: 4 mm	D1: 6 mm   D2: 5 mm
Description	Marker drill	Drilling extension	Trepine drill	
Material	Stainless steel			




- |   |   |  |  |
|---|---|--|--|
| 1<br><br>K5 - 9<br>driver<br>Ref. num.: B59   | 6<br><br>K5D - 15<br>driver<br>Ref. num.: DB515   | 11<br><br>K8N - 23<br>driver<br>Ref. num.: NB823 | 16<br><br>K1 - 15<br>driver<br>Ref. num.: B115           |
| 2<br><br>K5 - 15<br>driver<br>Ref. num.: B515 | 7<br><br>K8D - 23<br>driver<br>Ref. num.: DB823   | 12<br><br>K8N - 28<br>driver<br>Ref. num.: NB828 | 17<br><br>K9 - 23<br>adaptor<br>Ref. num.: B923          |
| 3<br><br>K8 - 23<br>driver<br>Ref. num.: B823 | 8<br><br>K8D - 28<br>adaptor<br>Ref. num.: DB828  | 13<br><br>K2 - 9<br>driver<br>Ref. num.: B29     | 18<br><br>K5 - 0<br>adaptor<br>Ref. num.: B50            |
| 4<br><br>K8 - 28<br>driver<br>Ref. num.: B828 | 9<br><br>K5N - 9<br>adaptor<br>Ref. num.: NB59    | 14<br><br>K2 - 15<br>driver<br>Ref. num.: B215   | 19<br><br>Parallel<br>pin<br>Ref. num.:<br>C8920 / C8927 |
| 5<br><br>K5D - 9<br>driver<br>Ref. num.: DB59 | 10<br><br>K5N - 15<br>adaptor<br>Ref. num.: NB515 | 15<br><br>K1 - 9<br>driver<br>Ref. num.: B19     | 20<br><br>R8T<br>ratchet<br>Ref. num.: C98               |



## SK-F

with  
Straight drills

Ref. num.: 10290



- |                                 |                                |                                  |                                  |                                  |                                  |                                  |                                  |                                  |
|---------------------------------|--------------------------------|----------------------------------|----------------------------------|----------------------------------|----------------------------------|----------------------------------|----------------------------------|----------------------------------|
| DR - 1.5<br>Ref. num.:<br>C9215 | DR - 2.0<br>Ref. num.:<br>C992 | DRF - 2.5<br>Ref. num.:<br>C9025 | DRF - 2.8<br>Ref. num.:<br>C9028 | DRF - 3.2<br>Ref. num.:<br>C9032 | DRF - 3.7<br>Ref. num.:<br>C9037 | DRF - 4.0<br>Ref. num.:<br>C9040 | DRF - 4.5<br>Ref. num.:<br>C9045 | DRF - 5.5<br>Ref. num.:<br>C9055 |
|---------------------------------|--------------------------------|----------------------------------|----------------------------------|----------------------------------|----------------------------------|----------------------------------|----------------------------------|----------------------------------|






- |   |   |  |  |
|---|---|--|--|
| 1<br><br>K5 - 9<br>driver<br>Ref. num.: B59   | 6<br><br>K5D - 15<br>driver<br>Ref. num.: DB515   | 11<br><br>K8N - 23<br>driver<br>Ref. num.: NB823 | 16<br><br>K1 - 15<br>driver<br>Ref. num.: B115           |
| 2<br><br>K5 - 15<br>driver<br>Ref. num.: B515 | 7<br><br>K8D - 23<br>driver<br>Ref. num.: DB823   | 12<br><br>K8N - 28<br>driver<br>Ref. num.: NB828 | 17<br><br>K9 - 23<br>adaptor<br>Ref. num.: B923          |
| 3<br><br>K8 - 23<br>driver<br>Ref. num.: B823 | 8<br><br>K8D - 28<br>adaptor<br>Ref. num.: DB828  | 13<br><br>K2 - 9<br>driver<br>Ref. num.: B29     | 18<br><br>K5 - 0<br>adaptor<br>Ref. num.: B50            |
| 4<br><br>K8 - 28<br>driver<br>Ref. num.: B828 | 9<br><br>K5N - 9<br>adaptor<br>Ref. num.: NB59    | 14<br><br>K2 - 15<br>driver<br>Ref. num.: B215   | 19<br><br>Parallel<br>pin<br>Ref. num.:<br>C8920 / C8927 |
| 5<br><br>K5D - 9<br>driver<br>Ref. num.: DB59 | 10<br><br>K5N - 15<br>adaptor<br>Ref. num.: NB515 | 15<br><br>K1 - 9<br>driver<br>Ref. num.: B19     | 20<br><br>R8T<br>ratchet<br>Ref. num.: C98               |

## SK-C



with  
Conical drills

Ref. num.: 10291






































- |                                 |                                |                                  |                                  |                                  |                                  |                                  |                                  |                                  |
|---------------------------------|--------------------------------|----------------------------------|----------------------------------|----------------------------------|----------------------------------|----------------------------------|----------------------------------|----------------------------------|
| DR - 1.5<br>Ref. num.:<br>C9215 | DR - 2.0<br>Ref. num.:<br>C992 | DRC - 2.5<br>Ref. num.:<br>C9125 | DRC - 2.8<br>Ref. num.:<br>C9128 | DRC - 3.2<br>Ref. num.:<br>C9132 | DRC - 3.7<br>Ref. num.:<br>C9137 | DRC - 4.0<br>Ref. num.:<br>C9140 | DRC - 4.5<br>Ref. num.:<br>C9145 | DRC - 5.5<br>Ref. num.:<br>C9155 |
|---------------------------------|--------------------------------|----------------------------------|----------------------------------|----------------------------------|----------------------------------|----------------------------------|----------------------------------|----------------------------------|










































 L6 Ø 2.0 Ref. num.: C10126	 L6 Ø 2.5 Ref. num.: C101256	 L6 Ø 2.8 Ref. num.: C101286	 L6 Ø 3.2 Ref. num.: C101326	 L6 Ø 3.7 Ref. num.: C101376	 L6 Ø 4.0 Ref. num.: C101406	 L6 Ø 4.5 Ref. num.: C101456
 L8 Ø 2.0 Ref. num.: C10128	 L8 Ø 2.5 Ref. num.: C101258	 L8 Ø 2.8 Ref. num.: C101288	 L8 Ø 3.2 Ref. num.: C101328	 L8 Ø 3.7 Ref. num.: C101378	 L8 Ø 4.0 Ref. num.: C101408	 L8 Ø 4.5 Ref. num.: C101458
 L10 Ø 2.0 Ref. num.: C101210	 L10 Ø 2.5 Ref. num.: C1012510	 L10 Ø 2.8 Ref. num.: C1012810	 L10 Ø 3.2 Ref. num.: C1013210	 L10 Ø 3.7 Ref. num.: C1013710	 L10 Ø 4.0 Ref. num.: C1014010	 L10 Ø 4.5 Ref. num.: C1014510
 L11,5 Ø 2.0 Ref. num.: C101211	 L11,5 Ø 2.5 Ref. num.: C1012511	 L11,5 Ø 2.8 Ref. num.: C1012811	 L11,5 Ø 3.2 Ref. num.: C1013211	 L11,5 Ø 3.7 Ref. num.: C1013711	 L11,5 Ø 4.0 Ref. num.: C1014011	 L11,5 Ø 4.5 Ref. num.: C1014511
 L13 Ø 2.0 Ref. num.: C101213	 L13 Ø 2.5 Ref. num.: C1012513	 L13 Ø 2.8 Ref. num.: C1012813	 L13 Ø 3.2 Ref. num.: C1013213	 L13 Ø 3.7 Ref. num.: C1013713	 L13 Ø 4.0 Ref. num.: C1014013	 L13 Ø 4.5 Ref. num.: C1014513

## SKD Stopper drills kit

Ref. num.: 104101





 L6 Ø 2.0 Ref. num.: C10126	 L6 Ø 2.5 Ref. num.: C105256	 L6 Ø 2.8 Ref. num.: C105286	 L6 Ø 3.2 Ref. num.: C105326	 L6 Ø 3.7 Ref. num.: C105376	 L6 Ø 4.0 Ref. num.: C105406	 L6 Ø 4.5 Ref. num.: C105456
 L8 Ø 2.0 Ref. num.: C10128	 L8 Ø 2.5 Ref. num.: C105258	 L8 Ø 2.8 Ref. num.: C105288	 L8 Ø 3.2 Ref. num.: C105328	 L8 Ø 3.7 Ref. num.: C105378	 L8 Ø 4.0 Ref. num.: C105408	 L8 Ø 4.5 Ref. num.: C105458
 L10 Ø 2.0 Ref. num.: C101210	 L10 Ø 2.5 Ref. num.: C1052510	 L10 Ø 2.8 Ref. num.: C1052810	 L10 Ø 3.2 Ref. num.: C1053210	 L10 Ø 3.7 Ref. num.: C1053710	 L10 Ø 4.0 Ref. num.: C1054010	 L10 Ø 4.5 Ref. num.: C1054510
 L11,5 Ø 2.0 Ref. num.: C101211	 L11,5 Ø 2.5 Ref. num.: C1052511	 L11,5 Ø 2.8 Ref. num.: C1052811	 L11,5 Ø 3.2 Ref. num.: C1053211	 L11,5 Ø 3.7 Ref. num.: C1053711	 L11,5 Ø 4.0 Ref. num.: C1054011	 L11,5 Ø 4.5 Ref. num.: C1054511
 L13 Ø 2.0 Ref. num.: C101213	 L13 Ø 2.5 Ref. num.: C1052513	 L13 Ø 2.8 Ref. num.: C1052813	 L13 Ø 3.2 Ref. num.: C1053213	 L13 Ø 3.7 Ref. num.: C1053713	 L13 Ø 4.0 Ref. num.: C1054013	 L13 Ø 4.5 Ref. num.: C1054513



**SKD-C**  
Conical stopper  
drills kit  
Ref. num.: 104105



3,75


## PLATFORM 3.75

# Prosthetic Accessories


One platform for P1 and P7

Internal hex: Perfect connection for abutments

Platform shifting




## Healing caps



Product code	H - 3.75, 2 mm	H - 3.75, 3 mm	H - 3.75, 4 mm	H - 3.75, 5 mm	H - 3.75, 6 mm
Ref. number	M10472	M10473	M10474	M10475	M10476
Dimensions	D: 4,7 mm H: 2 mm	D: 4,7 mm H: 3 mm	D: 4,7 mm H: 4 mm	D: 4,7 mm H: 5 mm	D: 4,7 mm H: 6 mm
Material	Titanium 6AL-4V				
Instructions	Suitable for all implant diameters. ⚠ Recommended tightening torque 15 Ncm for the screw.				

## Wide healing caps



[New!]



[New!]




[New!]




[New!]



Product code	HW - 3.75, 3 mm	HW - 3.75, 5 mm	HW7 - 3.75, 3 mm	HW7 - 3.75, 5 mm	HW8 - 3.75, 3 mm	HW8 - 3.75, 5 mm
Ref. number	M10583	M10585	M10703	M10705	M10803	M10805
Dimensions	D: 5,8 mm H: 3 mm	D: 5,8 mm H: 5 mm	D: 7,0 mm H: 3 mm	D: 7,0 mm H: 5 mm	D: 8,0 mm H: 3 mm	D: 8,0 mm H: 5 mm
Material	Titanium 6AL-4V					
Instructions	Suitable for all implant diameters.  Recommended tightening torque 25 Ncm for the screw.					

## Narrow healing caps



Product code	HN - 3.75, 3 mm	HN - 3.75, 5 mm
Ref. number	M10383	M10385
Dimensions	D: 3,8 mm H: 3 mm	D: 3,8 mm H: 5 mm
Material	Titanium 6AL-4V	
Instructions	Suitable for all implant diameters. ⚠ Recommended tightening torque 25 Ncm for the screw.	




## Open tray transfers



Product code	T1 - 3.75, 15 mm	T5 - 3.75, 11 mm	T4 - 3.75, 15 mm <i>slim</i>	T6 - 3.75, 11 mm <i>slim</i>	T8 - 3.75, 15 mm <i>slim</i>
Ref. number	M7615	M7611	M7615H	M7611S	M7615S
Dimensions	D: 4,5 mm H: 15 mm	D: 4,5 mm H: 11 mm	D: 3,6 mm H: 15 mm	D: 3,6 mm H: 11 mm	D: 3,6 mm H: 15 mm
Material	Stainless steel				
Instructions	⚠ Recommended tightening torque 25 Ncm for the screw.				


## Closed tray transfers



Product code	T1 - 3.75, 9 mm	T7 - 3.75, 9 mm <i>slim</i>	T2 - 3.75, 9 mm	T3 - 3.75, 9 mm
Ref. number	M779	M779S	M78	Transfer: M79 Plastic: M80
Dimensions	H: 9 mm	H: 9 mm	H: 9 mm	H: 9 mm
Material	Stainless steel			Titanium 6AL-4V/ Plastic
Instructions	⚠ Recommended tightening torque 25 Ncm for the screw.			



## Analog




Product code	A1 - 3.75
--------------	-----------

Ref. number	M08
-------------	-----

Dimensions	H: 12 mm
------------	----------


Material	Stainless steel
----------	-----------------

Instructions	 Suitable for all diameters of P1/P7
--------------	---

▶  
T5  
transfer


◀  
A1  
analog

## Straight abutments




Product code	S1 - 3.75, 5 mm	S1 - 3.75, 7 mm	S1 - 3.75, 9 mm	S1 - 3.75, 12 mm	S1 - 3.75, 15 mm
Ref. number	M115	M117	M119	M1112	M1115
Dimensions	H: 5 mm	H: 7 mm	H: 9 mm	H: 12 mm	H: 15 mm
Material	Titanium 6AL-4V				
Instructions	⚠ Recommended tightening torque 25 Ncm for the screw.				

## Straight narrow abutments




Product code	S1N - 3.75, 7 mm	S1N - 3.75, 9 mm
Ref. number	M11N7	M11N9
Dimensions	H: 7 mm	H: 9 mm
Material	Titanium 6AL-4V	
Instructions	⚠ Recommended tightening torque 25 Ncm for the screw.	

## Straight wide abutment





Product code	S1W - 3.75, 9 mm	S1W - 3.75, 12 mm
Ref. number	M11W9	M11W12
Dimensions	H: 9 mm	H: 12 mm
Material	Titanium 6AL-4V	
Instructions	⚠ Recommended tightening torque 25 Ncm for the screw.	

## Anatomic straight abutments



Product code	S1A - 3.75, 1 mm	S1A - 3.75, 2 mm	S1A - 3.75, 3 mm
Ref. number	M121	M122	M123
Dimensions	H: 1 mm	H: 2 mm	H: 3 mm
Material	Titanium 6AL-4V		
Instructions	⚠ Recommended tightening torque 25 Ncm for the screw.		

## Angular abutments 15°



Product code S2 - 3.75, 15° - 11 mm

S2L - 3.75, 15° - 15 mm

S2S - 3.75, 15° - 15 mm *slim*

Ref. number M131511

M131515

M141515

Dimensions H: 11 mm

H: 15 mm

H: 15 mm


Material

Titanium 6AL-4V

Instructions

⚠ Recommended tightening torque 25 - 30 Ncm for the screw.  
Two screw included in the package.

## Angular abutments 25°



Product code S2 - 3.75, 25° - 11 mm

S2L - 3.75, 25° - 15 mm

S2S - 3.75, 25° - 15 mm *slim*

Ref. number M132511

M132515

M142515

Dimensions H: 11 mm

H: 15 mm

H: 15 mm


Material

Titanium 6AL-4V

Instructions

⚠ Recommended tightening torque 25 - 30 Ncm for the screw.

## Angular abutments 35° - 45°



Product code S2 - 3.75, 35° - 11 mm

S2 - 3.75, 45° - 11 mm

Ref. number M133511

M134511

Dimensions H: 11 mm

H: 11 mm


Material

Titanium 6AL-4V

Instructions

⚠ Recommended tightening torque 25 - 30 Ncm for the screw.

## Anatomic angular abutments 15°




Product code	S2A - 3.75, 15° - 1 mm	S2A - 3.75, 15° - 2 mm	S2A - 3.75, 15° - 3 mm
Ref. number	M15151	M15152	M15153
Dimensions	H: 1 mm	H: 2 mm	H: 3 mm
Material	Titanium 6AL-4V		
Instructions	⚠ Recommended tightening torque 25 - 30 Ncm for the screw.		


## Anatomic angular abutments 25°



Product code	S2A - 3.75, 25° - 1 mm	S2A - 3.75, 25° - 2 mm	S2A - 3.75, 25° - 3 mm
Ref. number	M15251	M15252	M15253
Dimensions	H: 1 mm	H: 2 mm	H: 3 mm
Material	Titanium 6AL-4V		
Instructions	⚠ Recommended tightening torque 25 - 30 Ncm for the screw.		



## Anatomic zirconium abutments



Product code	S1AZ - 3.75, 1 mm	S1AZ - 3.75, 2 mm	S1AZ - 3.75, 3 mm
Ref. number	M161	M162	M163
Dimensions	H: 1 mm	H: 2 mm	H: 3 mm
Material	Zirconium		
Instructions	⚠ Recommended tightening torque 25 Ncm for the screw.		

## Anatomic angular zirconium abutments 15°



Product code	S2AZ - 3.75, 15° - 1 mm	S2AZ - 3.75, 15° - 2 mm	S2AZ - 3.75, 15° - 3 mm
Ref. number	M17151	M17152	M17153
Dimensions	H: 1 mm	H: 2 mm	H: 3 mm
Material	Zirconium		
Instructions	⚠ Recommended tightening torque 25 Ncm for the screw.		

## Anatomic angular zirconium abutments 25°



Product code	S2AZ - 3.75, 25° - 1 mm	S2AZ - 3.75, 25° - 2 mm	S2AZ - 3.75, 25° - 3 mm
Ref. number	M17251	M17252	M17253
Dimensions	H: 1 mm	H: 2 mm	H: 3 mm
Material	Zirconium		
Instructions	⚠ Recommended tightening torque 25 Ncm for the screw.		



## Anatomic zirconium abutments with titanium base



Product code	S1AZT - 3.75, 1 mm	S1AZT - 3.75, 2 mm	S1AZT - 3.75, 3 mm
Ref. number	M161T	M162T	M163T
Dimensions	H: 1 mm	H: 2 mm	H: 2 mm
Material	Zirconium / Titanium 6AL-4V		
Instructions	⚠ Recommended tightening torque 25 Ncm for the screw.		

## Anatomic angular zirconium abutments with titanium base 15°




Product code	S2AZT - 3.75, 15° - 1 mm	S2AZT - 3.75, 15° - 2 mm	S2AZT - 3.75, 15° - 3 mm
Ref. number	M17151T	M17152T	M17153T
Dimensions	H: 1 mm	H: 2 mm	H: 3 mm
Material	Zirconium / Titanium 6AL-4V		
Instructions	⚠ Recommended tightening torque 25 Ncm for the screw.		

## Anatomic angular zirconium abutments with titanium base 25°





Product code	S2AZT - 3.75, 25° - 1 mm	S2AZT - 3.75, 25° - 2 mm	S2AZT - 3.75, 25° - 3 mm
Ref. number	M17251T	M17252T	M17253T
Dimensions	H: 1 mm	H: 2 mm	H: 3 mm
Material	Zirconium / Titanium 6AL-4V		
Instructions	⚠ Recommended tightening torque 25 Ncm for the screw.		

## Anatomic straight peek abutments




Product code	S1A - Peek - 3.75, 1 mm	S1A - Peek - 3.75, 2 mm	S1A - Peek - 3.75, 3 mm
Ref. number	M181	M182	M183
Dimensions	H: 1 mm	H: 2 mm	H: 3 mm
Material	Peek-Classix polymer		
Instructions	<p>⚠ The peek abutments can be used for maximum 1 year. Recommended tightening torque 15 Ncm for the screw.</p>		

## Anatomic angular peek abutments




Product code	S2A - Peek - 3.75, 15° - 1 mm	S2A - Peek - 3.75, 15° - 2 mm	S2A - Peek - 3.75, 15° - 3 mm
Ref. number	M19151	M19152	M19153
Dimensions	H: 1 mm	H: 2 mm	H: 3 mm
Material	Peek-Classix polymer		
Instructions	<p>⚠ The peek abutments can be used for maximum 1 year. Recommended tightening torque 15 Ncm for the screw.</p>		

## Anatomic angular peek abutments




Product code	S2A - Peek - 3.75, 25° - 1 mm	S2A - Peek - 3.75, 25° - 2 mm	S2A - Peek - 3.75, 25° - 3 mm
Ref. number	M19251	M19252	M19253
Dimensions	H: 1 mm	H: 2 mm	H: 3 mm
Material	Peek-Classix polymer		
Instructions	<p>⚠ The peek abutments can be used for maximum 1 year. Recommended tightening torque 15 Ncm for the screw.</p>		

## Abutments for casting



Product code	S1PNH - 3.75, 11 mm	S1PN - 3.75, 11 mm
Ref. number	M82	M81
Dimensions	H: 11 mm	H: 11 mm
Material	Plastic	
Instructions	⚠ Recommended tightening torque 25 Ncm for the screw.	

## Abutment for casting



Product code	S1PT - 3.75, 11 mm	S1PCH - 3.75, 11 mm	S1PC - 3.75, 11 mm
Ref. number	M83	M85	M84
Dimensions	H: 11 mm	H: 11 mm	H: 11 mm
Material	Plastic / Titanium	Plastic / Chrome - Cobalt	Plastic / Chrome - Cobalt
Instructions	Melting range: > 900 °C	Melting range: > 1290 °C - 1380 °C	

## Abutments for bars





Product code	S1T - 3,75 / 0.5 mm	S1T - 3,75 / 1.5 mm	S1T - 3,75 / 2.5 mm	S1T - 3,75 / 3.5 mm
Ref. number	M8605	M8615	M8625	M8635
Dimensions	H: 0.5 mm	H: 1.5 mm	H: 2.5 mm	H: 3.5 mm
Material	Plastic / Titanium			
Instructions	⚠ Recommended tightening torque 25 Ncm for the screw.			

## Scanning abutment for 3D



Product code	SPC - 3.75, 11 mm
Ref. number	M107
Dimensions	11 mm
Material	Peek-Classix polymer

## Titanium base for the zirconium abutment




Product code	STC - 3.75 with hex	STCN - 3.75 without hex
Ref. number	M108	M109
Dimensions		
Material	Titanium 6AL-4V	
Instructions	Intended to be used for custom casting prosthetic restorations on single or multiple implants	Intended to be used for custom casting prosthetic restorations only on multiple implants

## Analog and screws



Product code	A1 - 3.75	S1a	S1aa
Ref. number	M08	110	110/b
Dimensions	H: 12 mm	H: 8 mm	H: 8 mm
Material	Stainless steel	Titanium 6AL-4V	
Instructions	Suitable for all diameters of P1/P7	For laboratory use	Standard abutment screw

## Individual block for milling




Product code	SIB - 3.75
Ref. number	M111
Dimensions	H: 25 mm
Material	Titanium 6AL-4V
Instructions	The titanium block abutment suitable for individual CAD CAM system

## Supported systems




*Dental Wings, 3shape, AmannGirrbach, Zirkonzahn and Exocad are trademarks or registered trademarks of their respective owners: Dental Wings LLC., 3Shape A/S, Amann Gribbach AG, Zirkonzahn USA Inc., Exocad GmbH.*

Ask our colleagues for detailed instructions and download links to library data.



## Ball attachment




Product code	S3 - 3.75, 1 mm	S3 - 3.75, 2 mm	S3 - 3.75, 3 mm	S3 - 3.75, 4 mm	S3 - 3.75, 5 mm	S3 - 3.75, 6 mm
Ref. number	M201	M202	M203	M204	M205	M206
Dimensions	1 mm	2 mm	3 mm	4 mm	5 mm	6 mm
Material	Titanium 6AL-4V					
Instructions	⚠ Recommended tightening torque 25 Ncm for the screw.					

## Angular ball attachment 17°



Product code	S3 - S7 - 3.75, 17° - 1 mm	S3 - S7 - 3.75, 17° - 2 mm	S3 - S7 - 3.75, 17° - 3 mm
Ref. number	M32171 / C4205	M32172 / C4205	M32173 / C4205
Dimensions	H: 1 mm	2 mm	3 mm
Material	Titanium 6AL-4V		
Instructions	⚠ Recommended tightening torque 25 Ncm for the screw. The basic package contains 0,5 mm ball attachment part.		



## Angular ball attachment 30°



Product code	S3 - S7 - 3.75, 30° - 1 mm	S3 - S7 - 3.75, 30° - 2 mm	S3 - S7 - 3.75, 30° - 3 mm
Ref. number	M32301 / C4205	M32302 / C4205	M323013 / C4205
Dimensions	1 mm	2 mm	3 mm
Material	Titanium 6AL-4V		
Instructions	⚠ Recommended tightening torque 25 Ncm for the screw. The basic package contains 0,5 mm ball attachment part.		



## S3 - S7 Ball Attachment parts



The package contains only the standrad silicone cap!

## Abutments for immediate loading



Product code	S4 - 3.75, 0.5 mm	S4 - 3.75, 1.5 mm	S4 - 3.75, 2.5 mm
Ref. number	M2605	M2615	M2625
Dimensions	H: 0.5 mm	H: 1.5 mm	H: 2.5 mm
Material	Titanium 6AL-4V		
Instructions	⚠ Recommended tightening torque 25 Ncm for the screw.		

## Angular abutments for immediate loading 17°



Product code	S4 - S7 - 3.75, 17° - 1 mm	S4 - S7 - 3.75, 17° - 2 mm	S4 - S7 - 3.75, 17° - 3 mm
Ref. number	M32171 / C43	M32172 / C43	M32173 / C43
Dimensions	H: 1 mm	H: 2 mm	H: 3 mm
Material	Titanium 6AL-4V		
Instructions	⚠ Recommended tightening torque 25 Ncm for the screw.		


## Angular abutments for immediate loading 30°



Product code	S4 - S7 - 3.75, 30° - 1 mm	S4 - S7 - 3.75, 30° - 2 mm	S4 - S7 - 3.75, 30° - 3 mm
Ref. number	M32301 / C43	M32302 / C43	M32303 / C43
Dimensions	H: 1 mm	H: 2 mm	H: 3 mm
Material	Titanium 6AL-4V		
Instructions	⚠ Recommended tightening torque 25 Ncm for the screw.		




[30°]




[0°]



[17°]



## S-Lock abutment - S35



Product code	S35 - 3.75, 1 mm	S35 - 3.75, 2 mm	S35 - 3.75, 3 mm	S35 - 3.75, 4 mm	S35 - 3.75, 5 mm	S35 - 3.75, 6 mm
Ref. number	M1121	M1122	M1123	M1124	M1125	M1126
Dimensions	H: 1 mm	H: 2 mm	H: 3 mm	H: 4 mm	H: 5 mm	H: 6 mm


Material

Titanium 6AL-4V

Instructions

⚠ Recommended tightening torque 30 Ncm for the screw.

## S-Lock angular abutment 17°



Product code	S35 - S7 - 3.75, 17° - 1 mm	S35 - S7 - 3.75, 17° - 2 mm	S35 - S7 - 3.75, 17° - 3 mm
Ref. number	M32171 / C113	M32172 / C113	M32173 / C113
Dimensions	H: 1 mm	H: 2 mm	H: 3 mm



Material

Titanium 6AL-4V

Instructions

⚠ Recommended tightening torque 30 Ncm for the screw.

## S-Lock angular abutment 30°



Product code	S35 - S7 - 3.75, 30° - 1 mm	S35 - S7 - 3.75, 30° - 2 mm	S35 - S7 - 3.75, 30° - 3 mm
Ref. number	M32301 / C113	M32302 / C113	M32303 / C113
Dimensions	H: 1 mm	H: 2 mm	H: 3 mm

Material

Titanium 6AL-4V

Instructions

⚠ Recommended tightening torque 30 Ncm for the screw.

## Accessories for S-Lock



**T1-S35 Transfer**  
Ref. num.:  
C114



**A1-S35 Analog**  
Ref. num.:  
C115



**K1 - 9 mm**  
Ref. num.:  
B19



**K1 - 15 mm**  
Ref. num.:  
B115



**K2 - 9 mm**  
Ref. num.:  
B29



**K2 - 15mm**  
Ref. num.:  
B215



**K21 Hand tool**  
Ref. num.:  
B21



**K9 - 23 mm**  
Ref. num.:  
B923

**Metal cap**  
Ref. num.:  
C116



**Laboratory**  
Ref. num.:  
C117



**Extra soft**  
Ref. num.:  
C118




**Soft**  
Ref. num.:  
C119



**Standard**  
Ref. num.:  
C120




**Strong**  
Ref. num.:  
C121



**Protective disk**  
Ref. num.:  
C122




The package contains all the caps illustrated above.



S35 recommended below 20 degrees, S35 - 17° / S35 - 30° recommended above 25 degrees of divergence between implant.

## Easy-Fix abutments - S5



Product code	S5 - 3.75, 1 mm	S5 - 3.75, 2 mm	S5 - 3.75, 3 mm	S5 - 3.75, 4 mm	S5 - 3.75, 5 mm	S5 - 3.75, 6 mm
Ref. number	M571	M572	M573	M574	M575	M576
Dimensions	1 mm	2 mm	3 mm	4 mm	5 mm	6 mm
Material	Titanium 6AL-4V					
Instructions	Recommended tightening torque 25 Ncm for the screw.					

## Easy-Fix angular abutments 17°



[New!]



[New!]



[New!]



Product code	S5 - S7 - 3.75, 17° - 1 mm	S5 - S7 - 3.75, 17° - 2 mm	S5 - S7 - 3.75, 17° - 3 mm
Ref. number	M32171 / C441	M32172 / C441	M32173 / C441
Dimensions	H: 1 mm	H: 2 mm	H: 3 mm
Material	Titanium 6AL-4V		
Instructions	Recommended tightening torque 25 Ncm for the screw. The basic package contains 1 mm S5 part.		

## Easy-Fix angular abutments 30°



[New!]



[New!]



[New!]



Product code	S5 - S7 - 3.75, 30° - 1 mm	S5 - S7 - 3.75, 30° - 2 mm	S5 - S7 - 3.75, 30° - 3 mm
Ref. number	M32301 / C441	M32302 / C441	M32303 / C441
Dimensions	H: 1 mm	H: 2 mm	H: 3 mm
Material	Titanium 6AL-4V		
Instructions	Recommended tightening torque 25 Ncm for the screw. The basic package contains 1 mm S5 part.		




## Accessories for Easy-Fix



T1-S5 Transfer  
Ref. num.: C58



A1-S5 Analog  
Ref. num.: C59



Metal cap  
Ref. num.: C50

Lab silicone  
Ref. num.: C54



Hard silicone  
Ref. num.: C51b


Medium silicone  
Ref. num.: C52

Soft silicone  
Ref. num.: C53

The package contains all the caps illustrated above.




K2 hand tool



S5 overdenture  
abutment



## S5 - S7 Easy-Fix parts




1 mm  
Ref. num.: C441

2 mm  
Ref. num.: C442

## Smart-Lock abutments - S8

H



Product code	S8 - 3.75, 1 mm	S8 - 3.75, 2 mm	S8 - 3.75, 3 mm	S8 - 3.75, 4 mm	S8 - 3.75, 5 mm	S8 - 3.75, 6 mm
Ref. number	M451	M452	M453	M454	M455	M456
Dimensions	H: 1 mm	H: 2 mm	H: 3 mm	H: 4 mm	H: 5 mm	H: 6 mm

Material

Titanium 6AL-4V

Instructions

⚠ Recommended tightening torque 25 Ncm for the screw.

## Smart-Lock angular abutment - S8 - 15°

15°



15°



Product code	S8 - 3.75, 15° - 1.5 mm	S8 - 3.75, 15° - 3 mm
Ref. number	M461515	M46153
Dimensions	H: 1.5 mm	H: 3 mm


Material

Titanium 6AL-4V

Instructions

⚠ Recommended tightening torque 25 Ncm for the screw.

## Package contains:



## Accessories for Smart-Lock abutment



T1-S8 Transfer  
Ref. num.:  
C48



A1-S8 Analog  
Ref. num.:  
C49




K4 Hand tool  
Ref. num.:  
C56

Insertion part




Magic tool  
Ref. num.:  
C55

Removal part




K4 hand tool




S8 overdenture  
abutment

## Smart-Lock overdenture abutment



## Straight multi-unit abutment - S6



Product code	S6 - 3.75, 0.5 mm	S6 - 3.75, 1.5 mm	S6 - 3.75, 2.5 mm
Ref. number	M3105	M3115	M3125
Dimensions	H: 0.5 mm	H: 1.5 mm	H: 2.5 mm
Material	Titanium 6AL-4V		
Instructions	⚠ Recommended tightening torque 25 Ncm for the screw.		

## Accessories for S6



◀ H - S6/S7  
Healing cap  
Ref. num.:  
C33



◀ T1 - S6/S7  
Transfer  
Ref. num.:  
C34



◀ A1 - S6/S7  
Analog  
Ref. num.:  
C35




◀ T - S6/S7  
Titanium sleeve  
Ref. num.:  
C36



◀ PH - S6/S7  
Plastic sleeve with hex  
Ref. num.:  
C37




◀ P - S6/S7  
Plastic sleeve without hex  
Ref. num.:  
C38



# Prosthetic options


A

Final restoration ▶



B

◀ Final restoration




## Multi-base angular abutment - S7 - 10°




Product code	S7 - 3.75, 10° - 1 mm	S7 - 3.75, 10° - 2 mm
Ref. number	M32101 / C40	M32102 / C40
Dimensions	1 mmH:	H: 2 mm
Material	Titanium 6AL-4V	
Instructions	⚠ Recommended tightening torque 25 Ncm for the screw.	

## Multi-base angular abutment - S7 - 17°



Product code	S7 - 3.75, 17° - 1 mm	S7 - 3.75, 17° - 2 mm	S7 - 3.75, 17° - 3 mm
Ref. number	M32171 / C40	M32172 / C40	M32173 / C40
Dimensions	H: 1 mm	H: 2 mm	H: 3 mm
Material	Titanium 6AL-4V		
Instructions	⚠ Recommended tightening torque 25 Ncm for the screw.		

## Multi-base angular abutment - S7 - 30°



Product code	S7 - 3.75, 30° - 1 mm	S7 - 3.75, 30° - 2 mm	S7 - 3.75, 30° - 3 mm
Ref. number	M32301 / C40	M32302 / C40	M32303 / C40
Dimensions	H: 1 mm	H: 2 mm	H: 3 mm
Material	Titanium 6AL-4V		
Instructions	⚠ Recommended tightening torque 25 Ncm for the screw.		

## Accessories for Multi-base abutments




◀ H - S6/S7  
Healing cap  
Ref. num.:  
C33



◀ T1 - S6/S7  
Transfer  
Ref. num.:  
C34



◀ A1 - S6/S7  
Analog  
Ref. num.:  
C35




◀ T - S6/S7  
Titanium sleeve  
Ref. num.:  
C36



◀ PH - S6/S7  
Plastic sleeve with hex  
Ref. num.:  
C37




◀ P - S6/S7  
Plastic sleeve without hex  
Ref. num.:  
C38



## S7 screw type abutment



◀ For more information  
check the video





## The-One multi-unit abutments - S16




Product code	S16 - 3.75, 1 mm	S16 - 3.75, 2 mm	S16 - 3.75, 3 mm	S16 - 3.75, 4 mm	S16 - 3.75, 5 mm
Ref. number	M641	M642	M643	M644	M645
Dimensions	H: 1 mm	H: 2 mm	H: 3 mm	H: 4 mm	H: 5 mm
Material	Titanium 6AL-4V				
Instructions	⚠ Recommended tightening torque 25 Ncm for the screw.				

## The-One angular multi-unit abutments - S17 - 17°




Product code	S17 - 3.75, 17° - 1 mm	S17 - 3.75, 17° - 2 mm	S17 - 3.75, 17° - 3 mm	S17 - 3.75, 17° - 4 mm	S17 - 3.75, 17° - 5 mm
Ref. number	M65171	M65172	M65173	M65174	M65175
Dimensions	H: 1 mm	H: 2 mm	H: 3 mm	H: 4 mm	H: 5 mm
Material	Titanium 6AL-4V				
Instructions	⚠ Recommended tightening torque 25 Ncm for the screw.				

## The-One angular multi-unit abutments - S17 - 30°



Product code	S17 - 3.75, 30° - 1 mm	S17 - 3.75, 30° - 2 mm	S17 - 3.75, 30° - 3 mm	S17 - 3.75, 30° - 4 mm	S17 - 3.75, 30° - 5 mm
Ref. number	M65301	M65302	M65303	M65304	M65305
Dimensions	H: 1 mm	H: 2 mm	H: 3 mm	H: 4 mm	H: 5 mm
Material	Titanium 6AL-4V				
Instructions	⚠ Recommended tightening torque 25 Ncm for the screw.				

## The-One angular multi-unit abutments - S17 - 45°



Product code	S17 - 3.75, 45° - 1 mm	S17 - 3.75, 45° - 2 mm	S17 - 3.75, 45° - 3 mm	S17 - 3.75, 45° - 4 mm	S17 - 3.75, 45° - 5 mm
Ref. number	M65451	M65452	M65453	M65454	M65455
Dimensions	H: 1 mm	H: 2 mm	H: 3 mm	H: 4 mm	H: 5 mm
Material	Titanium 6AL-4V				
Instructions	⚠ Recommended tightening torque 25 Ncm for the screw.				

## Accessories for S16 / S17



T1 - S16/S17  
Transfer  
Ref. num.: C67



P - S16/S17  
Plastic sleeve  
without hex  
Ref. num.: C71



H - S16/S17  
Healing cap  
Ref. num.: C66




A1 - S16/S17  
Analog  
Ref. num.: C68




T - S16/S17  
Titanium sleeve  
Ref. num.: C69



PC - S16/S17  
Chrome Cobalt sleeve  
Ref. num.: C70





All solutions for a **perfect smile.**

Drill  
Surgical  
Prosthetic  
Implant  
Quality  
Skills  
Warranty  
About  
Surgical kits  
Titanium  
Zirconium  
Knowledge  
Favorable prices  
Surgical instruments  
Ratchet Hex Driver  
Bone type  
SGS De  
Wide range of sizes  
Anatomic  
Plastic sleeve  
Education  
Long  
Peek  
Surgical tools  
Ball attachment  
Angular  
Crilling  
T  
Conus platform  
ernal Hex  
g cap  
s6  
9s  
Drill



**SGS Dental Implant System Holding, reg No: FL-00023**

FL-9494, Schaan, Landstrasse 27  
Liechtenstein  
Tel.: 00423 233 5050, 00423 233 5051  
Fax: 00423 233 5052

---

**SGS International Kft./Ltd. European Logistic Center**

H-1047 Budapest, Károlyi István u. 1-3.  
Tel: +36 1 328 0427

---

**[www.sgs-dental.com](http://www.sgs-dental.com)**